CATALOG #3300 Catalogue # 3300

**GEN-PROBE®**  
**PACE®**  
Use only with the GEN-PROBE PACE Systems and AMPLIFIED Chlamydia Trachomatis Assay

**CONTENTS:**  
50 SPECIMEN COLLECTION KITS  
Each kit contains:  
Two control swabs  
One tube containing 0.9 mL GEN-PROBE transport medium

Store at Room Temperature  
For *in vitro* diagnostic use

**GEN-PROBE®**  
**PACE®**  
Use only with the GEN-PROBE PACE Systems and AMPLIFIED Chlamydia Trachomatis Assay

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**GEN-PROBE®**  
**PACE®**  
Use only with the GEN-PROBE PACE Systems and AMPLIFIED Chlamydia Trachomatis Assay

**SPECIMEN  
COLLECTION  
KITS FOR  
ENDOCERVICAL SPECIMENS**

Use only with the GEN-PROBE PACE Systems and AMPLIFIED Chlamydia Trachomatis Assay

**SPECIMEN  
COLLECTION  
KITS FOR  
ENDOCERVICAL SPECIMENS**

**For *in vitro* diagnostic use**  
Store at Room Temperature

CATALOG #3300

**Specimen Collection and Handling of  
Endocervical Specimens for  
C. trachomatis (PACE and AMPLIFIED) and/or  
N. gonorrhoeae (PACE only) Testing**

1. Remove excess mucus from cervical os and surrounding mucosa using one of the swabs provided. Discard this swab.
2. Insert second swab from collection kit 1 to 1.5cm into endocervical canal.
3. Rotate swab clockwise for 10 to 30 seconds in endocervical canal to ensure adequate sampling.
4. Withdraw swab carefully; avoid any contact with vaginal mucosa.
5. Prepare swab for transport.

**Preparation for Transport**

1. **Fully insert one swab** into the GEN-PROBE transport tube.
2. Snap off shaft at score line or cut shaft to fit tube. Use care to avoid splashing of contents.
3. Cap tube tightly.
4. For best results, transport tube to the laboratory at 2° to 25°C. Store at 2° to 25°C and test within 7 days of collection.
5. If not tested within 7 days:
  - a. PACE: Vortex specimen, remove swab, and store at -20° to -70°C.
  - b. AMPLIFIED: Store at -20° to -70°C.**Do not remove swab.**

**Note: Specimens collected with this system cannot be used for culture. Only swabs with the GEN-PROBE specimen collection system should be used for specimen collection.**

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